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4D Laser Mask

【First use must do】

- First try on the face for 1 minute, no tingling, redness, discomfort and then use normally.

【Core usage frequency】

- It can be used daily, up to 2 times, and the cumulative use time should not exceed 30 minutes per day.

*For other detailed information, please refer to the instruction manual.

【Safety Taboos】(Must-read)

- Do not use if wearing medical electronic devices such as pacemakers, electrocardiogram monitors, artificial heart-lung machines, or infusion pumps.
- Do not use in the bathroom or near water sources. Keep away from containers filled with water.
- Keep away from heating sources such as stoves or air conditioning vents.
- Do not use if the device, remote control, or power cord is damaged.
- Do not disassemble or modify the device by yourself.

【Quick steps】


Before use, the skin must be cleaned and kept dry.


Put on the mask and secure it on your face to ensure a comfortable fit.


On the controller, press and hold the power button for 2 seconds to turn on the machine. The lights will start to illuminate. For the first use, test for 1 minute. If there is no discomfort, you can use it normally.

Press the timer button briefly to adjust the care time. (The default setting upon startup is T1:15 minutes. You can adjust for different time periods: T1:15 minutes, T2:20 minutes, T3:30 minutes)

There are 3 modes. Press "+" to increase, "-" to decrease. Adjust to the mode that suits you. Complete the care according to the preset time and the machine will automatically shut down after the completion.

After use: Remove the mask and proceed with your daily skincare routine.