



SLIMAX ΣΥΣΚΕΥΗ ΑΝΤΙΓΗΡΑΝΤΙΚΗΣ ΦΡΟΝΤΙΔΑΣ

【Φόρτιση και χρήση】

- Είσοδος: 5V 1A – χρησιμοποιήστε τον παρεχόμενο φορτιστή
- Βεβαιωθείτε ότι η επαφή φόρτισης είναι σωστά συνδεδεμένη
- Η συσκευή δεν είναι αδιάβροχη – μην τη βυθίζετε σε νερό

【Σημαντικές αντενδείξεις / προειδοποιήσεις】

- Μην χρησιμοποιείτε σε διαβρωτικά περιβάλλοντα
- Μην αποσυναρμολογείτε τη συσκευή. Να τη χειρίζεστε με προσοχή.
- Αν το δέρμα παρουσιάσει ενόχληση, μειώστε την ένταση ή διακόψτε άμεσα τη χρήση

【Γρήγορα βήματα】

- 01 Καθαρίστε το δέρμα σας. Τοποθετήστε το φύλλο βαμβακιού στην κεφαλή της συσκευής, ασφαλίστε το με τον δακτύλιο στερέωσης και προσθέστε επαρκή ποσότητα από το serum ή το καθαριστικό σας.
- 02 Ενεργοποιήστε τη συσκευή και επιλέξτε από τις παρακάτω λειτουργίες περιποίησης:

EPT Facial Care – Daily Care
(Καθημερινή περιποίηση προσώπου)**EPT Facial Care – Deep Treatment**
(Βαθιά θεραπεία προσώπου)**SNE Body Care –**
Περιποίηση σώματος όπως για ώμους
και αυχένα**MAS Head Care –**
Περιποίηση κεφαλιού

- 03 Ρυθμίστε την ένταση. Για πρώτη χρήση, συνιστάται επίπεδο 1–2.
- 04 Εκτελέστε απαλή περιποίηση σε πρόσωπο / σώμα / κεφάλι:
Για το πρόσωπο: ανοδικές κινήσεις (lifting massage) όπως φαίνεται στην εικόνα
Για το κεφάλι: τεχνική μασάζ σε σχήμα "Z"
- 05 Μετά τη χρήση, απενεργοποιήστε τη συσκευή και τοποθετήστε την στη θήκη αποθήκευσης.
Χτυπήστε ελαφρά το δέρμα για να απορροφηθεί το serum.

JOVS



QUICK START GUIDE CARD

Slimax Anti-aging Device

【Charging and Use】

- Input: 5V 1A, use the included power supply for charging
- Ensure that the charging contact is in good contact before charging
- Not waterproof. Do not immerse in water

【Important Taboos】

- Do not use in corrosive environments
- Do not disassemble or drop the device by yourself
- If the skin feels uncomfortable, reduce the gear or stop using it immediately.

【Quick steps】

- 01 Prepare your skin for cleaning and use it with the fixed ring, cotton pad, and collagen essence.
02 Turn on the device and select the care mode:

EPT Facial Care - Daily Care



EPT Facial Care - Deep Treatment



SNE Body Care - Such as Shoulder and Neck Care



MAS Head Care



- 03 Adjust the intensity of the gears. For the first use, it is recommended to set to 1-2 gears.
04 Perform gentle facial/ body/ head care: such as the upward pulling massage technique for the face daily, or using the "Z" shape massage technique on the head.
05 After use, turn off the device, and put it back into the storage box. And pat the skin to allow the essence to be absorbed.