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Neck mask

【First use must do】

- Light Sensitivity Test: Apply it to the forearm and expose it to light for 10 minutes. If there is no redness, itching or other reactions, you can use it normally.

【Important Precautions】(Note)

- Pregnant women, people with photosensitivity, and patients with skin diseases are prohibited from using
- Those wearing medical devices such as pacemakers are prohibited from using
- Do not use in the bathroom or near water sources as it is not waterproof
- Do not use while charging
- Do not attempt to disassemble or modify it yourself

【Quick steps】

Clean the skin: Clean the neck and dry it. Keep it makeup-free and dry.

Assemble: Attach the main unit and the controller, adjust the straps and then wear it.

Power On: Press the power button to turn on. The default setting is low gear; short press to switch from low to medium to high mode.

Automatic Shutdown: Once turned on, it starts a 20-minute timer. When the timer is over, it automatically shuts down.

Post-Routine Care: Remove the instrument and continue with your daily skincare routine; it is recommended to use 3-5 times per week.